


BTL-08 SD3 DIAG 30.32 PM 4.75.0000
Aut

16/06/2025 11:16:04


10mm/s
20mm/mV
50Hz ADS 1.5s




HR = 63



HR = 63



HR = 63



CE 0187

CE 0187

UVJ	P+	P-	Q	R	S	R	S	T+	T-	ST	STs [uV/s]	SF	RR	63 l/min
161	-	-	-44	1089	-163	-	-	251	-	20	763	P	952 ms	
133	-	-	-	1276	-111	-	-	251	-	34	535	QRS	94 ms	
-	-	-	-	43	-86	385	-	17	-20	10	-223	PQ	144 ms	
-	-	-	-	136	-	-	-	-	-246	-26	-642	QT	378 ms	
VR	82	-139	-151	473	-170	-	-	123	-	0	460	QTc	387 ms	
VL	53	-	-51	751	-	-	-	130	-	22	136	axis P	15°	
VF	44	-	-	309	-952	-	-	439	-53	8	-217	axis QRS	45°	
146	-	-	-	1385	-942	-	-	198	-	72	1134	axis T	24°	
99	-	-	-	1206	-268	-	-	198	-	18	577			
83	-	-	-	1165	-146	-	-	190	-	18	392			
78	-	-	-	1134	-121	-	-	195	-	9	641			
62	-	-	-	1039	-61	-	-	174	-	21	386			

ECG DANS DES LIMITES NORMALES
RYTHME SINUSAL

A F Curly